



A safe, effective and convenient countermeasure against COVID-19

EcoProven Anti-Viral Mouth Wash / Spray is a safe, effective and convenient preventative measure against SARS-CoV-2 virus - that causes COVID-19. The EcoProven technology turns nature's most common solvent (water) into a powerful virucide. The formula is a proprietary blend of plant-based and all natural ingredients (ionic and anionic compounds) that laboratory testing has confirmed is effective at penetrating the fatty lipid structure that surrounds and protects the SARS-CoV-2 virus.

EcoProven helps prevent initial infection of the SARS-Cov-2 virus in the upper respiratory tract. As a spray application or mouth rinse, EcoProven covers the inside surface of the mouth and throat and acts quickly to neutralize the virus before it can embed itself and migrate deeper into the lungs.

EcoProven Anti-Viral Mouth Wash / Spray is also helpful in mitigating early stage infections and reducing symptoms. For individuals who have tested positive, the product helps to disrupt the on-set of COVID -19 symptoms. For individuals who are actively presenting symptoms, the product helps to stall the progression of the illness.

The potent biological effect of EcoProven can stay active for hours - sustained by the presence of moisture (water) in the upper respiratory tract. In this way, EcoProven as a mouth wash or spray is a very effective defence against COVID-19 and helps reduce the spread of the virus to others.

Physical Properties


| | | | |
|-----------------|------------------|--------------------|--------------------------------|
| Physical State: | Liquid at 70 F | pH: | 10.4 NEAT (in concentrate) |
| Color: | Clear | Foaming: | Extremely low foaming |
| Odor: | Odorless | Viscosity: | @20c > 225 - 375 (290 avg) cps |
| Taste: | Mint or Cinnamon | Pounds per gallon: | 8.52 |
| Flashpoint: | None at 212 F | Specific Gravity: | 1.029 |
| Freezing Point: | 32 F | Effective Period: | As long as water is present |

EPA Company ID#: **97213** EPA Facility ID #: **2935261** FDA Product #: **pending**

How EcoProven Works

EcoProven's technology combines with the moisture in your mouth and upper respiratory tract to become a hyperactive colloidal water solution. This solution contains countless microscopic ionic particles which discharge electrical charges in all directions. The virus body is "caught in the crossfire" of these ionic charges and the protective outer lipid membrane of the virus is shredded. Any virus that is in contact with EcoProven's powerful colloidal action is destroyed.

Once this outer protective layer is destroyed, the inner RNA is exposed and the combination of the colloidal action and antimicrobial surfactants quickly destroys the virus' RNA core as well.



Technology Approvals

- FDA Reg. # - pending
- CDC Reg. # - pending
- EPA Water Sense Certified
- CSA Sustainability Certified
- UL Green Guard Certified
- Green Seal Certified
- Approved & Used by United States Navy
- Approved & Used by United States Postal Service
- Approved by United States
Department of Agriculture

Protect from temperatures below 32 F and above 90 F
No special handling / storage procedures are needed

